


# QR-scanning



Open the rear facing camera on your mobile. On some devices, you need to tap the AI icon that looks like this: 

In older devices, you may need a QR reader app. Point the camera at the QR code. The device should detect the QR code and display a message.



Tap the notification to open the link associated with the QR code. Follow the on-screen instructions to complete the login.

## Link the code manually

A screenshot of a web form titled 'Register new kit'. It has a red header bar. Below the header is a text input field labeled 'Enter code' with the value 'abc123' entered. At the bottom of the form is a dark grey button labeled 'Register'.

Alternatively, you can visit [dev-selfsampling.hpvelimination.se/myprofile](https://dev-selfsampling.hpvelimination.se/myprofile) and connect your code manually.

## Successful linkage

A screenshot of a web page titled 'Register your kit'. It has a green header bar. Below the header is a green message box that says 'Your kit (kit code: abc123) is now linked to your identity. You can now proceed with the instructions and send the sample to us'. At the bottom of the page is a red button labeled 'Proceed to consent'.

Once your sample has been linked to your identity, you will see this message. Now you can take the test and send it to us!